**READ BEFORE USING PROCESSED  
HUMAN ALLOGRAFT TISSUE**

**THIS TISSUE WAS RECOVERED FROM A DECEASED  
DONOR WHOSE LEGAL NEXT OF KIN HAS GIVEN  
PERMISSION FOR THE BONE AND CONNECTIVE TISSUE TO  
BE DONATED OR FROM A FEMORAL HEAD FROM A  
CONSENTING LIVING DONOR. RECOVERY WAS  
PERFORMED IN ACCORDANCE WITH EUROPEAN  
DIRECTIVES.**

**PROCESSING AND PACKAGING WERE PERFORMED  
UNDER PHARMACEUTICAL CONDITIONS.**

**DESCRIPTION AND INDICATION FOR USE**

The Tissue establishment provides freeze-dried bone tissue for clinical use, in a variety of standard sizes to meet the needs of qualified health care professionals (e.g., physicians, dentists, and/or podiatrists) and impregnated with vancomycin (0.1 g vancomycin added per 1 cc bone). Processed human bone is used in a variety of surgical applications and in combination with prosthetic devices. The choice of type and size/dosage of allograft depends on the surgeon's preference and on the kind of defect.

**PACKAGING**

Freeze-dried tissues are aseptic impregnated with vancomycin and packaged in a bottle in a pouch, which is heat-sealed. The content of the vial is sterile. The outer pouch is labeled with barcode and human readable for the donor identification sequence of the SEC, additionally lot number, product code, expiration date, tissue description, size, quantity (volume) is also indicated on the pouch label.

Additional information is listed on the package label.

**CAUTION**

1. Every package is for single patient use only; do not use portions of the same allograft from a single package on multiple patients;
2. Preparation of the bone graft bed is important for allograft incorporation and bone formation; there must be no infections at the site of grafting;
3. Caution should be taken when the patient is allergic to antibiotics, because the allograft is impregnated with antibiotic (vancomycin).

**PRECAUTION**

Transmission of infectious diseases and onset of bacterial infections at the site of grafting may occur despite the extensive medical procedures and serological testing performed to select all tissue donors.

Any adverse outcomes attributable to this allograft must be reported promptly to the Tissue establishment.

**Donor Screening and Testing**

Before recovery, the donor's medical and social history was obtained to exclude medical conditions or disease processes that are contra indicated for tissue donation, in accordance with current EU directives. All tissue bank policies and procedures for donor screening, serological and microbiological testing, meet or exceed current standard guidelines. Contraindications for tissue donation include: presence of infectious disease, malignant disease, neurological degenerative diseases, diseases of unknown etiology, and exposure to toxic substances.

According to current guidelines, the Tissue Bank evaluates the infectious disease testing results, together with the consent documentation, donor medical history, risk behaviors, physical assessment, available medical records including previous medical history, laboratory test results existing autopsy or coroner reports, as well as information from other sources or records which may pertain to donor eligibility along with procurement test results. These data are sufficient to establish if the donor meets current eligibility criteria at the time of recovery, which are incompliance with current regulations.

Donor blood samples were tested at least for:

- HIV-1 /2 antibody
- Syphilis
- Hepatitis B surface antigen
- Hepatitis B core antibody
- Hepatitis C antibody
- HTLV (when applicable)

Tissue cultures, collected at recovery, were tested for aerobic, anaerobic microbial contaminants and yeasts.

Name and address of the testing laboratories, list and interpretation of all required serological and microbiological tests, list of the documents reviewed as part of the relevant medical records and name of the person(s) or establishment determining the eligibility of this human tissue, are on file by Tissue Bank. Those records are available upon request of the competent authority and retained for 35 years.

This allograft tissue has been determined to be suitable for transplantation.

**STERILITY**

Discard unused allograft portions or in excess, in accordance with recognized procedures for medical waste disposal.

This allograft is guaranteed sterile in an unopened and undamaged package.

Do not use this allograft if:

- packaging is damaged or tampered,
- the package label or the barcode is damaged, unreadable or missing,
- bottle containing the allograft is broken,
- the allograft has expired.

If any of the above conditions exist or are suspected, please notify the Tissue establishment.

**Sterilization with Gamma Irradiation**

The tissue is processed in facilities with classified environments. The tissue is terminally sterilized and virus inactivated by gamma irradiation.

**FREEZE-DRIED BONE AND DEMINERALIZED FREEZE-DRIED BONE**

Lyophilization (freeze-drying) lowers the total water content to <5% in freeze-dried allograft.

**STORAGE**

Store freeze-dried allograft at room temperature (15°C-25°C).

It is the responsibility of the transplant facility or clinicians to keep the tissue for clinical use in the recommended storage conditions before transplant.

**RECONSTITUTION/REHYDRATION PROCEDURE**

To obtain the best clinical results and prevent graft failure, the procedure and recommendations listed below should be followed.

**Preparation for Use**

All reconstitution techniques should be performed using aseptic/sterile technique.

It is common surgical practice to reconstitute tissue in a normal sterile saline solution.

Patient sensitivity to antibiotics used to reconstitute allograft tissues should be checked prior to use.

Use one new solution for each allograft. Solvent should be sufficient to cover the tissue completely.

**Freeze-dried tissue should be reconstituted for at least one hour. Before use, the allograft should be carefully soaked in normal sterile saline solution.** If tissue must be shaped, drilled or cut, it is recommended to reconstitute it for a longer time, so to increase flexibility depending on the surgeon's needs.

If tissue is not used directly after reconstitution, it must be stored at 2°C - 8°C in sealed containers until its use. If the reconstituted tissue is not used within 2 hours after reconstitution, it must be discarded.

Do not reconstitute the allograft in a solution or medium warmer than 52 °C.

**Use standard aseptic/sterile technique to open the package.**

Please, open freeze-dried tissue package as follows:

- Peel off the flaps of the pouch;
- Open the vial;
- Put the allograft into a sterile basin for reconstitution;
- Soak completely the allograft into reconstitution solution.\*See above

**References:**

- European Directive 2004/23/EC - 31 March 2004
- European Directive 2006/17/EC - 08 February 2006
- European Directive 2006/86/EC - 24 October 2006

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