

METAPIN 3mm - 5mm METAPIN HOLDER



Titanium membrane tacks

Product Features

- Unbreakable and unbendable
- Sterile and ready to use
- Single packaged
- Tack placement instrument with precision guide



METAPIN HOLDER

METAPIN 3 mm

METAPIN 5 mm



METAPIN 3mm-5mm METAPIN HOLDER



The membrane tacks METAPIN are made out of Titanium in order to guarantee high resistance and maximum biocompatibility. METAPIN are specifically recommended for the placement of both resorbable and non-resorbable membranes even in areas of the jaw where the bone density is higher.

METAPIN HOLDER is an autoclavable bone tack placement instrument proposed with a protective cap that must be removed during the phase of pickup, and can be used in order to complete the stabilization process by dabbing on the head of the positioned tack, when required. METAPIN HOLDER features an exclusive precision guide that guarantees an improved level of stability of the tack during the phase of insertion and stabilization, thus allowing a more precise placement and a stronger grip.

METPIN are available in two different lengths: 3 and 5 mm, in a ready-to-use sterile package of 3 tacks each. The sterile packaging improves the safety of the clinical management reducing the time connected to the process of sterilization and cleaning.

Ref. 4710 METAPIN HOLDER 1pc

Ref 4720 METAPIN 3 mm. 3 pcs, Sterile.

Ref. 5030 METAPIN 5 mm. 3 pcs, Sterile.

OFICINA EN ESPAÑA

Balmes 190, 4º
08006 Barcelona
T +34 932 380 735

salugraft@salugraftdental.com
www.salugraftdental.com

OFICINA EN PORTUGAL

Loures Business Park, En 115, Lote 10
Sul Escritório Nº 1
2660-515 São Julião do Tojal
T +351 961 637 375

