


INTEROSS[®]

ANORGANIC CANCELLOUS BONE GRAFT GRANULES



WHAT IS INTEROSS[®]

InterOss[®] is a natural hydroxyapatite bone grafting material for bone regenerative application in dentistry. Made from a proven multi-step purification process, InterOss[®] is an osteoconductive material for bone regeneration. InterOss[®] is available in granule form. It is supplied sterile and is dedicated for single use.



InterOss[®] Vial Type

Size	0.25-1.0 mm			
Reference	IOSG025	IOSG050	IOSG100	IOSG200
Weight	0.25g	0.50g	1.00g	2.00g
Volume	0.54cc	1.08cc	2.16cc	4.32cc

Size	1-2 mm		
Reference	IOLG050	IOLG100	IOLG200
Weight	0.25g	0.50g	1.00g
Volume	0.54cc	1.08cc	2.16cc



SIGMAGRAFT


THE BONE GRAFT SUBSTITUTE SPECIALISTS

INDICATIONS FOR USE


InterOss[®] granules are recommended for:

- ✔ Augmentation or reconstructive treatment of the alveolar ridge
- ✔ Filling of periodontal defects
- ✔ Filling of defects after root resection, apicoectomy, and cysectomy
- ✔ Filling of extraction sockets to enhance preservation of the alveolar ridge
- ✔ Elevation of the maxillary sinus floor
- ✔ Filling of periodontal defects in conjunction with products intended for Guided Tissue Regeneration (GTR) & Guided Bone Regeneration (GBR)
- ✔ Filling of peri-implant defects in conjunction with products intended for GBR


BENEFITS



Micropore (<10µm)



Mesopore (10-100µm)



Granules

- ✔ Biocompatible and osteoconductive making it ideal for bone regeneration uses
- ✔ Hydrophilic properties allow rapid absorption of patient's blood making it easy to apply
- ✔ Interconnected network of macro and micro pores providing an ideal environment for cell attachment and bone growth

Contact Us

Address: SigmaGraft, Inc. 575 Sally Place Fullerton, CA 92831, USA
Phone: (714) 525-0114, Toll Free: (888) 499-0114
E-Mail: info@sigmagraft.com