



# ZMATRIX™

*porcine peritoneum collagen membrane*

## PERFECTLY SOFT

*Sometimes you need a membrane pliable enough to contour to its environment but substantial enough to maintain its architecture. Introducing Zmatrix™ - a natural, porcine peritoneum collagen membrane with a perfectly soft consistency that drapes without the usual self-adherence experienced with other natural collagen membranes.*

- ▶ *Designed to drape without adhering to itself*
- ▶ *Natural peritoneum collagen structure provides for elasticity*
- ▶ *Processed to preserve extracellular components including laminin, fibronectin, elastin, and glycosaminoglycans\**





# ZMATRIX™

porcine peritoneum collagen membrane

Zmatrix™ is a natural, native collagen membrane; cross-linking chemicals and agents are unnecessary. Proprietary processing technology allows preservation of collagen as well as extracellular components including laminin, fibronectin, elastin, and glycosaminoglycans.\*

## Zmatrix™ Membrane Sizes

Sold Individually



15 x 20 mm

20 x 30 mm

30 x 40 mm

Part No. ZM1520

Part No. ZM2030

Part No. ZM040

\*Hoganson DM, Owens GE, O'Doherty EM, Bowley CM, Goldman SM, Harilal DO, Neville CM, Kronengold RT, Vacanti JP. Preserved extracellular matrix components and retained biological activity in decellularized porcine mesothelium. *Biomaterials*. 2010, 27: 6934-6940.