



OraGRAFT™

Unicortical Block

Clinical Overview The OraGraft portfolio is comprised of allograft implants generally used for sinus augmentation, craniofacial reconstruction and correction of periodontal and ridge defects.

Applications Alveolar Ridge Augmentation

-
- Features & Benefits**
- **Safety:** Sterile with a Sterility Assurance Level (SAL) of 10^{-6} .
 - **Convenient Storage:** Ambient temperature with 5 year shelf life stored with LifeNet Health's proprietary Preservon® technology.
 - **Ease of Use:** Can be affixed to surgical site using standard bone screws.
 - **Optimal Remodeling:** Cortico/Cancellous blocks allow for space maintenance and remodeling due to trabecular architecture.
 - **Match to Defect:** Easy to trim, can be shaped using standard instrumentation.
-





OraGRAFT Unicortical Block

Ambient Temperature Storage (10°C - 37°C), 5 year shelf life

Size	Order Code
10 x 15 x 8 mm	UB101508
10 x 25 x 8 mm	UB102508

Instructions for use available at [LifeNetHealth.org/IFU](https://www.lifenethealth.org/IFU)